



| | |
|---------------|-------------------|
| Product Code | 45186 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279451863 |
| PARAF | 977817016 |

Description

FLAVIA LIFT ARMCHAIRS WITH 2 MOTORS

2 motors lift armchair with standard or extra lifting footrest: it can reach the bed position providing a relaxing experience. Its backrest and footrest can operate separately to any position. Extra lifting footrest allows 8 cm lifting adjustment higher than flat position.

- Color: brown
- 2 motors
- Fireproof: NO
- Footrest: extra

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position.

Power supply: 100-240 V - 50/60 Hz.

Multilingual manual: GB, FR, IT, ES.

Technical Specifications

- Seat
Width: 47 cm
Depth: 50 cm
Height: 47 cm
- Backrest height: 74 cm
- Chair
Width: 78 cm
Length (reclined): 171.5 cm
- Load: 150 Kg (331 lbs)