

STOOL with backrest and ring - black



Product Code	45220
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279452204
---------	---------------

Description

Fire proof chair with padding and covering made of simulated leather.
Height adjustable. Plastic base with 5 castors and chromed ring.

Height adjustable 60-86 cm
Maximum load: 110 kg

Made in Italy.