



|               |                        |
|---------------|------------------------|
| Product Code  | 24526                  |
| Unit of sale  | box of 40 pcs          |
| Minimum order | 1                      |
| Type          | IVD - Professional Use |
| Class         |                        |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2498736       |
| CND        | W0105090102   |
| EAN/UPC    | 6936983144704 |

Description

HELICOBACTER PYLORI TEST - box of 40 tests - PROFESSIONAL

Rapid test for qualitative detection of Helicobacter pylori antibodies in human whole blood, serum or plasma.  
Easy to perform and interpret.  
Quick results in 10 minutes.

It includes 40 test cassettes, 40 droppers, 1 buffer.

Multilingual box and manual: GB, FR, IT, ES, DE, HU.

Relative sensitivity: 96.8%; relative specificity: 93.0% and accuracy: 94.6%.

WARNING: SHORT EXPIRY DATE FOR TESTS. Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.