


CRP TEST - semi quantitative



Product Code	24527
Unit of sale	box of 20 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	1808212
CND	W0103020607
EAN/UPC	8023279245271
PARAF	922289487

Description

CRP TEST - box of 20 tests

Measurement of CRP in whole blood, plasma or serum has two purposes:

- 1) hs-CRP levels help in assess the risk of cardiovascular and peripheral vascular disease. Higher hs-CRP values have highest risk of cardiovascular disease.
 - 2) CRP levels help identify and keep track of infections and diseases that cause inflammation (lymphoma, lupus, giant cell arteritis, osteomyelitis, rheumatoid arthritis,).
- use fingerstick whole blood or serum/plasma
 - results within 7 minutes
 - room temperature storage
 - built in performance control and reference band

It includes 20 cassette tests with droppers. For professional use.

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.