



Product Code	45290
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279452907

Description

MULTIFUNCTIONAL WALL UNITS - STANDARD SIZE

Colour: white

Dimensions: 103 x 20 x h 40 cm

Multi-purpose wall unit to organize dispensers of disposable products.
All components and slots are modular and can be organized to match professional user needs.
Powder coated galvanized metal pull-up door with vertical opening.
Supplied with automatic soap dispenser, cups dispenser and 5 slots.

Standard accessories

Dispenser for disposable cups 1
Automatic soap dispenser 1
Slot for gloves or wipes 2
Slot for masks 1
Roll holder 2