


MULTI-PURPOSE WALL UNIT - large - white



Product Code	45295
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279452952

Description

MULTIFUNCTIONAL WALL UNITS - LARGE SIZE

Colour: white
Dimensions: 153 x 20 x h 40 cm

Multi-purpose wall unit to organize dispensers of disposable products.
All components and slots are modular and can be organized to match professional user needs.
Powder coated galvanized metal pull-up door with vertical opening.
Supplied with automatic soap dispenser, cups dispenser and 7 slots.

Standard accessories

- Dispenser for disposable cups 1
- Automatic soap dispenser 1
- Slot for gloves or wipes 3
- Slot for masks 1
- Slot for napkins or cannula 1
- Roll holder 2