



Product Code	24531
Unit of sale	box of 25 pcs
Minimum order	1
Type	IVD - Professional Use
Class	
UK-REP	Yes
RDM (NSIS)	2450597
CND	W0102039013
EAN/UPC	8023279245318

Description

FOB RAPID TEST - box of 25 tests

Accurate test to qualitatively detect FOB levels as low Bas 40 ng/ml. It reads results in 5 minutes.
It allows early patient treatment decisions.
Simple procedure and interpretation.

It includes 25 cassettes and specimen collection tubes with extraction buffer, multilingual manual and box: GB, FR, IT, ES, PT, DE, GR, RU.

WARNING: SHORT EXPIRY DATE FOR TESTS. Most tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Technical Specifications

Sensitivity: 97.6%
Specificity: 99.3%
Accuracy: 98.9%