

JACKET - cotton/polyester - unisex L light blue



Product Code	21454
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	8023279214543
---------	---------------

Description

JACKET - UNISEX - 60% COTTON - 40% POLYESTER

- Size: L
- Colour: light blue
- UNISEX
- Short sleeves
- 2 front and 1 breast pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.
Stock in 5 colours, other 22 colours available in 60 days.
White jacket 200 g/m2, coloured jacket 185 g/m2.

Technical Specifications

Composition: 60% cotton - 40% polyester
Weight: 200 g/m2 (white), 185 g/m2 (coloured)
Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2
Construction: warp: 40 ± 3 weft: 19 ± 3
Breaking strength - warp: 100 kg ±3%
Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)
Fastness to washing: 4
Treatments and finishing: crease resist
Dimensional stability - washing at 90°C: 3%