

# MOBILE UNIT GE4LP under wash-basin 49 cm - any colour



Colori disponibili  
Available colours  
Couleurs disponibles



Product Code	45453
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279454536

## Description

### MOBILE UNITS

Model: under washbasin

Available in 45 days, please specify colour with order

Mobile space is a programme of safe and ergonomic technical furniture for medical, physiotherapy and professional beauty use. These units feature wheels and offer the same levels of functionality as trolleys. The visible edges are radiating and covered in shockproof material. The steel handles provide good grip without dangerous edges. Stability is optimised by the inner and outer sides, the extra-thick backs and invisible wheel positioning. The steel runners of the drawers with rocker arm closing ensure smooth opening while preventing undesired opening during transport. The independent worktop makes it possible to fully pull out the units for easier floor cleaning and position the trolleys according to need. It is available with or without washbasin (left or right). The broad choice of unit configurations and colour range make it possible to personalise interiors according to individual requirements. Covered in shock-proof material. Easy removal by front pull handle. Invisible wheels.

Made in Italy.