




Product Code	24547
Unit of sale	box of 40 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	1752484
CND	W0105090401
EAN/UPC	8023279245479

Description

MALARIA RAPID TEST - box of 40 tests

A rapid chromatographic immunoassay test for the qualitative detection of circulating plasmodium falciparum in whole blood. Compared with traditional thick or thin blood smears microscopic analysis, the correlation between the two systems is over 99,0%.

Package includes test devices, buffer, disposable specimen dropper and multilingual manual (GB, FR, IT, ES), lancets are not included.

For professional use.

WARNING: SHORT EXPIRY DATE FOR TESTS.
Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.