

ADJUSTABLE FOLDABLE ARM - without curtain



Product Code	45501
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279455014

Description

Swing arm curtain - Adjustable foldable arm (without curtains).

Swing arm curtain is adjustable horizontally (180°) and can be bent down to 90° when not in use.

Telescopic tubular arm extension: open 210 cm, closed 100 cm.

Rings allow curtains to slide smoothly.