

TROLLEY for 1 curtain - foldable - without curtain



Product Code	45505
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279455052
---------	---------------

Description

Trolley curtain - 1 curtain (without curtains).

Trolley made of stainless steel AISI 304 with 4 revolving castors (2 lockable).

It allows the use of 1 curtain 225 x h 145 cm.

Arm can be bent down to 90° when not in use and curtain remains tied up.