

TROLLEY for 2 curtains - foldable - without curtain



Product Code	45506
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279455069

Description

Trolley curtain - 2 curtains (without curtains).

Trolley made of stainless steel AISI 304 with 4 revolving castors (2 lockable).
It allows the use of 2 curtains (1 long 225 x h 145 cm and 1 short 175 x h 145 cm).
Arm can be bent down to 90° when not in use and curtain remains tied up.