

# CLOSED RINGS for use with wall arms and trolleys



|               |                   |
|---------------|-------------------|
| Product Code  | 45509             |
| Unit of sale  | box of 18 pcs     |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279455090 |

## Description

Set of 18 replacement rings for Trevira® curtain system

- Closed rings for wall arms and trolleys