

TREVIRA CURTAINS for arm - 225 x h180 cm - light blue



Product Code	45510
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279455106
PARAF	951864115

Description

Trevira® curtains

- Curtain for wall arms and cabin
- Size: 225 x h 180 cm
- Colour: light blue

Trevira® curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C.

Made in Italy.