

# TREVIRA CURTAINS for trolley - long 225 x h145 cm - light blue



|               |                   |
|---------------|-------------------|
| Product Code  | 45520             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279455205 |

## Description

Trevira® curtains

- Curtain for trolleys and windoes (Arm/Rail, Long trolley)
- Size: 225 x h 145 cm
- Colour: light blue

Trevira® curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C.

Made in Italy.