

# TREVIRA CURTAINS for trolley - short 175 x h145 cm - light blue



Product Code 45530

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279455304

## Description

Trevira® curtains

- Curtain for trolleys and windows (Arm/Rail, Short trolley)
- Size: 175 x h 145 cm
- Colour: light blue

Trevira® curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C.

Made in Italy.