

# TREVIRA CURTAINS for trolley - short 175 x h145 cm - blue



Product Code	45531
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279455311
---------	---------------

## Description

Trevira® curtains

- Curtain for trolleys and windoes (Arm/Rail, Short trolley)
- Size: 175 x h 145 cm
- Colour: blue

Trevira® curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C.

Made in Italy.