

4 WING SCREEN with castors - without curtains



Product Code	45563
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279455632
PARAF	965140041

Description

SCREEN WITH TREVIRA CURTAINS (not included)
Aluminium alloy tubular frame mounted on castors (Ø 50 mm).

- 4 wing screen with castors
- Size: 200xh 167.5 cm

Packing size: 176x57x13 cm

Made in Italy.