

TREVIRA CURTAINS for wing screen - blue



Product Code	45571
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279455717
PARAF	965140054

Description

CURTAINS IN TREVIRA 45xh 129 cm

Trevira curtains contain a permanent silver-based antimicrobial agent which is not affected by frequent washing.

UV resistant and flame-retardant material.

Disinfectable and washable at 60°.

Available in 6 colours.

Custom-made curtains available on request.

Multilingual manual: GB, FR, IT.

Made in Italy.

Made in Italy.