

3 WINGS SCREEN - white cloth



Product Code	45580
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2225411
CND	V0899
EAN/UPC	8055136921127

Description

Exclusive 3 wings ward screens with sturdy structure in oval tube (30 x 15 mm), painted with epoxidic powder metal grey.

- Width: 60 + 60 + 60 cm (3 sections)
- Height: 180 cm
- White washable cotton cloth curtains
- Equipped with 4 adjustable feet that ensure increased stability

Made in Italy.