

MODULAR TROLLEY painted steel with 1 drawer + 1 shelf



Product Code	45603
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279456035
GTIN	08023279456035

Description

MODULAR TROLLEYS SYSTEM

- Trolley with 1 small drawer, 1 compartment with shelf
- Painted steel

Modular trolleys made of painted steel or stainless steel, one- piece structure, with wide technopolymer plastic worktop with borders on 3 sides to be cleaned easily, no sharp corners, ergonomic tubular push-handle. Different drawers (small-yellow, medium-blue, large-red), depending on model, with epoxy painted panels and integrated handle grip, silent sliding with end-run block. Possibility to remove the drawers from the guides for an easy cleaning and disinfection. Central lock key. Overall bumpers and 4 antistatic Ø 125 mm castors, 2 with brakes.

Size: 67x63x104 cm.

Weight: 53 kg.

Made in Italy.