



Product Code	45656
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1418766
CND	V0880
EAN/UPC	8023279456561
GMDN	36801

## Description

### PEDIATRIC KART

Pediatric compact kart with 2 medium and 2 large drawers.

Size: 654 x 605 x h 975 mm

High quality aluminium multifunctional cart with superior features.

- work-top in high performance copolymer, with perimetric edge for liquid containment
- extruded aluminium frame, with superficial silver antioxidant treatment offering more lightness, resistance and reliability
- ergonomic stainless steel handles, fitted on both sides for easy handling
- drawer structure built in steel sheets, epoxy powder coated
- full pull-out sliding drawers
- interchangeable drawer position
- aluminium drawer handle with "hook/unhook" system to prevent accidental opening of drawers
- high drawers capacity
- centralized locking system with front key
- 6 integrated side pockets (3+3) to hold bottles or small items
- swivelling wheels Ø 125 mm (two with brake) in thermoplastic rubber with threadguards, hub on precision ball bearings.

## Technical Specifications

2 medium and 2 large drawers

Drawer size:

- medium h2= 2/10

Front width/height: 490xh 137 mm

Inner sizes: 432x448xh 131 mm

- large h3= 3/10

Front width/height: 490xh 210 mm

Inner sizes: 432x448xh 204 mm