



Product Code	45667
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	964315
CND	V0880
EAN/UPC	8023279456677
GMDN	36801

## Description

EMERGENCY MEDIUM KART - 5 drawers - red

Size: 804 x 616 x h 975 mm

High quality aluminium emergency cart with superior features.

- work-top in high performance copolymer, with perimetric edge for liquid containment
- extruded aluminium frame, with superficial silver antioxidant treatment offering more lightness, resistance and reliability
- ergonomic stainless steel handles, fitted on both sides for easy handling
- drawer structure built in steel sheets, epoxy powder coated
- full pull-out sliding drawers
- interchangeable drawer position
- aluminium drawer handle with "hook/unhook" system to prevent accidental opening of drawers
- high drawers capacity
- centralized locking system with single use seals
- 6 integrated side pockets (3+3) to hold bottles or small items
- swivelling wheels Ø 125 mm (two with brake) in thermoplastic rubber with threadguards, hub on precision ball bearings

## Technical Specifications

5 drawers: 2 small, 1 medium, 2 large

Drawer size:

- small h1= 1/10

Front width/height: 640xh 64 mm

Inner sizes: 582x448xh 58 mm

- medium h2= 2/10

Front width/height: 640xh 137 mm

Inner sizes: 582x448xh 131 mm

- large h3= 3/10

Front width/height: 640xh 210 mm

Inner sizes: 582x448xh 204 mm