



| | |
|---------------|----------------|
| Product Code | 45668 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 964315 |
| CND | V0880 |
| EAN/UPC | 8023279456684 |
| GMDN | 36801 |

Description

EMERGENCY MEDIUM KART - 4 drawers - orange

Size: 804 x 616 x h 975 mm

High quality aluminium emergency cart with superior features.

- work-top in high performance copolymer, with perimetric edge for liquid containment
- extruded aluminium frame, with superficial silver antioxidant treatment offering more lightness, resistance and reliability
- ergonomic stainless steel handles, fitted on both sides for easy handling
- drawer structure built in steel sheets, epoxy powder coated
- full pull-out sliding drawers
- interchangeable drawer position
- aluminium drawer handle with "hook/unhook" system to prevent accidental opening of drawers
- high drawers capacity
- centralized locking system with single use seals
- 6 integrated side pockets (3+3) to hold bottles or small items
- swivelling wheels Ø 125 mm (two with brake) in thermoplastic rubber with threadguards, hub on precision ball bearings

Technical Specifications

4 drawers: 2 medium, 2 large

Drawer size:

- medium h2= 2/10

Front width/height: 640xh 137 mm

Inner sizes: 582x448xh 131 mm

- large h3= 3/10

Front width/height: 640xh 210 mm

Inner sizes: 582x448xh 204 mm