



| | |
|---------------|----------------|
| Product Code | 45671 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 1418767 |
| CND | V0880 |
| EAN/UPC | 8023279456714 |
| GMDN | 36801 |

Description

VIEWER KART

Compact kart with slim LED viewer, 2 small drawers and 2 folder holders.

Size: 654 x 605 x h 975 mm

High quality aluminium multifunctional cart with superior features.

- work-top in high performance copolymer, with perimetric edge for liquid containment
- extruded aluminium frame, with superficial silver antioxidant treatment offering more lightness, resistance and reliability
- ergonomic stainless steel handles, fitted on both sides for easy handling
- drawer structure built in steel sheets, epoxy powder coated
- full pull-out sliding drawers
- interchangeable drawer position
- aluminium drawer handle with "hook/unhook" system to prevent accidental opening of drawers
- high drawers capacity
- centralized locking system with front key
- 6 integrated side pockets (3+3) to hold bottles or small items
- swivelling wheels Ø 125 mm (two with brake) in thermoplastic rubber with threadguards, hub on precision ball bearings.

Technical Specifications

2 small drawers and 2 folder holders

Drawer size:

- small h1= 1/10

Front width/height: 490xh 64 mm

Inner sizes: 432x448xh 58 mm

- folder holders h4= 4/10

Largeur frontale/hauteur: 490xh 283 mm

Dimensions internes: 432x448xh 270 mm