



Product Code	45672
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1418767
CND	V0880
EAN/UPC	8023279456721
GMDN	36801

Description

VIEWER KART 2

Compact kart with slim LED viewer, 2 small and 1 medium drawers, 10 compartments for folders (44 x 48 x h 43.8 cm).

Size: 654 x 605 x h 975 mm

High quality aluminium multifunctional cart with superior features.

- work-top in high performance copolymer, with perimetric edge for liquid containment
- extruded aluminium frame, with superficial silver antioxidant treatment offering more lightness, resistance and reliability
- ergonomic stainless steel handles, fitted on both sides for easy handling
- drawer structure built in steel sheets, epoxy powder coated
- full pull-out sliding drawers
- interchangeable drawer position
- aluminium drawer handle with "hook/unhook" system to prevent accidental opening of drawers
- high drawers capacity
- centralized locking system with front key
- 6 integrated side pockets (3+3) to hold bottles or small items
- swivelling wheels Ø 125 mm (two with brake) in thermoplastic rubber with threadguards, hub on precision ball bearings.

Technical Specifications

2 small and 1 medium drawers

Drawer size:

- small h1= 1/10

Front width/height: 490xh 64 mm

Inner sizes: 432x448xh 58 mm

- medium h2= 2/10

Front width/height: 490xh 137 mm

Inner sizes: 432x448xh 131 mm