



Product Code	45684
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279456844

Description

DRESSING TROLLEY - small

Dressing trolley is supplied with upper dispenser (10 tilting bins).

- antimicrobial epoxy painted body with key lock
- top surface with 4-sided rim
- top surface with built-in compartment
- ball-bearing full-extension sliding drawer with self-closing system
- built-in handle
- 4 swivelling castors Ø 125 mm - 2 with brake
- lateral waste container 15 l with lid and catheter holder
- size: small 760x714xh 1,022 mm

Made in Italy.