



Product Code	45711
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279457117

Description

MULTIFUNCTIONAL TROLLEY 65x88xh 100 cm

Solid and robust construction with steel sheet painted with antimicrobial epoxy powders and technopolymers for maximum lightness, maneuverability and ease of use.
 Upper part: large support surface in PST with ergonomic and integrated push handle. Push handle has a central mechanism for locking and unlocking wheels' brakes.
 Front part: retractable worktop support with "push and pull" opening.
 Four modular drawers with soft self-closure, ball bearings and central locking with key.
 Drawers are made of painted steel.

Left side: recessed compartment and retractable worktop made of technopolymer plastic; 10 L waste bin with knee opening and technopolymer catheter holder.

Right side: recessed compartment with three tilting pockets made of technopolymer material.
 Supplied with four Ø 100 mm twinned castors with automatic locking system.

Made in Italy.