



Product Code	45715
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279457155

Description

NEO PLUS CLINICAL TROLLEY

Clinical trolley made of aluminium, steel and ABS structure.

Upper part is equipped with stainless steel 304 three-sided rail, which prevents items from falling.

It is equipped with four castors Ø 100 mm (two with brake), one drawer with silent guide rails and ABS inserts, which can be freely organized, and card holder.

Size: 62.5x47.5x92 cm.

Max load: table top 80 kg, middle tray 15 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, PL,RO, CZ, HU, SE, FI, LT, HR, GR, BG.

Table board covered by transparent glass

Table top internal size: 57.5x43x1 cm

Card holder

Drawer inner size: 43x33.5x11 cm

Internal tray size: 47.5x36x5 cm

Internal bottom tray size: 56x42x2.5 cm

Standard accessories

- 2 ABS waste bins 7 l with lid
- 2 shelves
- 1 drawer with dividers
- 1 fixed tray with inner hollow depression