

OXYGEN CONCENTRATOR 10 L



Product Code	34572
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2582202
CND	Z12159004
EAN/UPC	6933257924451

Description

SILENT OXYGEN CONCENTRATOR 10 I

Medical grade, silent, small oxygen concentrators for home use.
Provide high oxygen concentration through nasal cannula.
Work continuously.

Supplied with filter, humidifier bottle, nasal cannula and manual GB, IT.
Available on gimaitaly.com: ES, SE, GR, TR, Arabic.

Technical Specifications

Flow range: 1 ~ 10 l/min
Oxygen concentration: 87 - 96 %
Output pressure: 150 kPa
Input power: 700 VA
Noise level (front - whole unit): 49 dBA
Size - weight: 41x37.5x66.8 cm - 24 kg
Packing size: 51x44x75 cm
Operating voltage: 230 V - 50 Hz
Optimal oxygen concentration: after 15 min.
Classification: Type BF, Class II equipment