



Product Code	45764
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2429909
CND	Z129085
EAN/UPC	8053736713180
UMDNS	11726
GMDN	41893

Description

Plate adapter for monitors trolley code 45760

- Type: D

Other brands adapter

- Brand: Colin
- Model: Press-Mate 8800

- Brand: Fukuda Denshi
- Model: Dynascope DS-8007 and DSA-81

- Brand: General Electric
- Model:
GE Dash 2000, 3000, 4000 monitors
Carescape B650, B450
Carescape B20, B40, B105, B125

- Brand: Masimo
- Model:
Root M series
Rad-97 M series

- Brand: Nellcor
- Model:
N-595, N-550 Pulse Oximeter
Bedside SpO2 PM100N

- Brand: Nihon Kohden
- Model:
Life Scope i, Life Scope L, BSM 2301, BSM 2303
Vismo PVM 2700, Vismo PVM 4000
Life Scope VS, BSM 3500, BSM 3700, Vismo PVM 2700
Life Scope P, BSM 4103, BSM 4113
Life Scope A, BSM 5105, BSM 5135
Life Scope TR, BSM 6301K, BSM 6501K, BSM 6701K
Life Scope G5, Life Scope N, OPV 1500
GF-210R Multigas Unit
CSM-1501, CSM-1502, BSM-6XXX, WEP-5200, GF-200 series

- Brand: Oridion
- Model: Capnostream™ 20 and Capnostream™ 20p

- Brand: Philips
- Model:
Efficia CM10-12-100-120-150
EarlyVue VS30

- Brand: Sanuwave
- Model: MIST and UltraMIST series

- Brand: Sonosite
- Model: S and M series

Made in Italy.