



Product Code	45792
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 8023279457926

Description

UTILITY TROLLEY 2 shelves + 1 drawer

Practical trolleys with stainless steel frame and two or three high quality ABS shelves.
All shelves provided with three-sided stainless steel siderails to prevent items from falling down.
Two push handles for easy moving.
Wear resistant, high quality silent castors, Ø 100 mm, two with brake.

Size: 73x48x90 cm
Drawer size: 50x42x12 cm
Max load: 60 kg (30 kg per shelf).