



| | |
|---------------|-------------------|
| Product Code | 45805 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279458053 |

Description

TER PLUS TROLLEY

With removable trays
 Stainless steel frame, mounted on castors Ø 80 mm.
 Size: 49x30xh 78 cm.
 Weight: 9 kg.
 Maximum load (x shelf): 10 kg.
 Delivered in kit form.

Made in Italy.