

# LANCART TROLLEY 60x40xh80 cm with drawer



Product Code	45808
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 8023279458084

## Description

Lancart trolley - small

- Size: 60x40xh 80
- Shelves: 2
- Drawer: yes

Stainless steel trolleys mounted on castors Ø 80 mm.  
Maximum trolley load: 60 kg.  
Range of optional accessories to customize your trolley.  
Drawer size: 29x30xh 7,5 cm  
Delivered in kit form.

Made in Italy.