

LANCART TROLLEY 70x50xh80 cm with drawer



Product Code 45809

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

EAN/UPC 8023279458091

Description

Lancart trolley - medium

- Size: 72x50xh 80
- Shelves: 2
- Drawer: yes

Stainless steel trolleys mounted on castors Ø 80 mm.

Maximum trolley load: 60 kg.

Range of optional accessories to customize your trolley.

Drawer size: 29x30xh 7,5 cm

Delivered in kit form.

Made in Italy.