



Product Code	45810
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279458107
PARAF	951840622

Description

Lancart trolley - small
 - Size: 60x40xh 80
 - Shelves: 2
 - Drawer: -

Stainless steel trolleys mounted on castors Ø 80 mm.
 Maximum trolley load: 60 kg.
 Range of optional accessories to customize your trolley.
 Delivered in kit form.

Made in Italy.