

LANCART TER TROLLEY 70x50xh80 cm - 3 shelves



Product Code	45813
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279458138

Description

Lancart Ter trolley - medium
- Size: 72x50xh 80
- Shelves: 3
- Drawer: -

Stainless steel trolleys mounted on castors Ø 80 mm.
Maximum trolley load: 60 kg.
Range of optional accessories to customize your trolley.
Delivered in kit form.

Made in Italy.