

OXYGEN CONCENTRATOR 5 L



Product Code	34582
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2969580
CND	Z12159004
EAN/UPC	6945040108245
PARAF	950975906

Description

Small and portable oxygen concentrator composed of air compressor, absorption tower of molecular sieve, sensor and flowrate controller.
The machine adopts PSA (pressure swing absorbers) method to directly concentrate medical oxygen with high concentrator from the air. Lightweight, easy to use and safe.
Stable performance and low noise level.
Wide application in clinic, hospital, oxygen bar, gym center...

Supplied with humidifier bottle, power cord, nasal cannula and user manual (GB, FR*, IT, ES*, DE*) * only on www.gimaitaly.com

Size: 508x260x h 530 mm.
Weight: 21 kg

Technical Specifications

Max. flowrate: 5l/min
Flow range: 0.5~5l/min
Oxygen concentration: 93% ±3%
Output pressure: 20kPa~50kPa
Input power: 450VA
Noise level: =55dB(A)
Operating voltage: AC220V ± 10%, 50Hz ±1 Hz
Type of protection: Class II, Type BF, IP21