

HYDRAULIC MAYO TABLE - 5 SPOKES



Product Code 45834

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP No

EAN/UPC 8023279458343

Description

HYDRAULIC MAYO TABLE with pump

Serving trolley used as instrument holder in operating rooms.
Not removable tray made of AISI 304 stainless steel.
Equipped with hydraulic pump for height adjustment.
Base made of aluminium with 5 spokes and swiveling castors Ø 50 mm.
Delivered in kit form.

Made in Italy.

Technical Specifications

S/S tray size: 70x45xh 1.8 cm
Adjustable height: 85 to 120 cm
Load: 40 kg
Weight: 14 kg
Packing size: 60x17xh 77 cm