



Product Code	45835
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
--------	----

EAN/UPC	8023279458350
GTIN	08023279458350

Description

HYDRAULIC MAYO TABLE with pump and removable tray

Height adjustable instruments trolley with removable tray in AISI 304 stainless steel, scotch brite finish.

Tray can be rotated through 360° with tightening knob.

Base in AISI 304 stainless steel tubular frame fixed to the hydraulic pump that allows height adjustment by foot pedal.

Three castors Ø 75 mm, one antistatic and two with brake.

Delivered in kit form.

Made in Italy.

Technical Specifications

S/S tray size: 70x45xh 1.8 cm

Adjustable height: 90 to 126 cm

Load: 20 kg

Weight: 17 kg

Packing size: 98x68xh 15 cm