



Product Code	45835
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
--------	----

EAN/UPC	8023279458350
---------	---------------

Description

HYDRAULIC MAYO TABLE with pump and removable tray

Height adjustable instruments trolley with removable tray in AISI 304 stainless steel, scotch brite finish.

Tray can be rotated through 360° with tightening knob.

Base in AISI 304 stainless steel tubular frame fixed to the hydraulic pump that allows height adjustment by foot pedal.

Three castors Ø 75 mm, one antistatic and two with brake.

Delivered in kit form.

Made in Italy.

Technical Specifications

S/S tray size: 70x45xh 1.8 cm

Adjustable height: 90 to 126 cm

Load: 20 kg

Weight: 17 kg

Packing size: 98x68xh 15 cm