



Product Code	45881
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2805588
CND	Z129099
EAN/UPC	8023279458817

Description

DOUBLE BOWL STAND - Stainless Steel

Aisi 202 stainless steel stand and base with 5 plastic wheels, one with brake.

Adjustable height from 60 to 95 cm.

Aisi 304 stainless steel bowl Ø 32 cm, 5 l capacity.

Delivered in kit form.

Base Ø 65 cm

Made in Italy.