

SPIRIT PORTABLE OXYGEN CONCENTRATOR 1 I



Product Code	34589
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2567515
CND	Z12159004
EAN/UPC	6933257924925
UMDNS	12873
PARAF	950975920

Description

SPIRIT PORTABLE OXYGEN CONCENTRATOR 1 I

Portable and silent oxygen concentrator designed to provide a reliable and efficient source of oxygen. It provides high concentration of oxygen to patients with respiratory disorders. Its compact design and ease of transport is ideal for those who require continuous oxygen therapy anywhere.

- 4 settings for regulating the amount of oxygen delivered
- detachable molecular sieve
- 0.02 s response
- oxygen generation by impulses
- modular lithium battery
- automatic warning system

Supplied with rechargeable lithium-ion battery, AC power supply + DC power cable, 1 set of felt filter (2 pieces), nasal cannula, carrying bag with protective case and adjustable shoulder strap, manual in GB, IT.

Technical Specifications

Flow range: 3.5 l/min
Oxygen concentration: 96~90%
Output pressure: 150 kPa
Operation noise level: = 55 dBA
Input power: 120 W
Voltage: AC 100-240 V, 50/60 Hz
Battery: DC 14.4 V
Power Supply: DC 12V-24 V
Optimal oxygen concentration: after 30 min
Size: 22x8xh 22 cm
Packing size: 32x32xh 21 cm
Weight: 2.2 kg (with battery)
Classification: Class II equipment, type BF applied part, IP20