

# SOILED LINEN TROLLEY - stainless steel AISI 202



Product Code	45890
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Yes
--------	-----

CH-REP	No
--------	----

EAN/UPC	8023279458909
---------	---------------

## Description

SOILED LINEN TROLLEY - stainless steel AISI 202

Linen trolley made of stainless steel, intended for collection and transport of soiled linen at hospital facility.  
Supplied with 4 castors and removable bag.

Size: 65x55xh 90 cm