

DOUBLE SOILED LINEN TROLLEY - stainless steel AISI 202



Product Code	45891
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Yes
CH-REP	No
EAN/UPC	8023279458916

Description

DOUBLE SOILED LINEN TROLLEY - stainless steel AISI 202

Linen trolley made of stainless steel, intended for collection and transport of soiled linen at hospital facility. Supplied with 4 castors and 2 removable bags.

Size: 94x42xh 90 cm