



Product Code	24590
Unit of sale	box of 20 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2037460
CND	W0105099099
EAN/UPC	8023279245905
PARAF	983782398

## Description

SALIVA COVID TEST - professional - box of 20

The COVID-19 Antigen Rapid Test (Oral Fluid) is a qualitative rapid chromatographic immunoassay for the qualitative detection of SARSCoV-2 Nucleocapsid (N) protein antigens in oral fluid specimens from individuals with suspected infection. It's a single-use test kit intended to detect the coronavirus SARS-CoV-2 through human oral fluid. This test is designed for professional use.

Research shows that the SARS-CoV-2 virus is stable in oral fluid for long periods of time at warm temperatures.

- fast results: 15 minutes
  - high accuracy
  - convenient procedure
  - easy interpretation
  - simple operation
  - relative specificity: 99,3%
- Storage temperature: 2-30°C

Shelf life: 24 months (Gima stock may vary from 6 to 22 months)

## Standard accessories

Test devices  
Collection devices (funnel, tube and tube tip)  
Buffer  
Package insert  
Biosafety bag  
Procedure card  
Multilingual box: GB, FR, IT, ES, DE, PT.