



Product Code	45900
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279459005

## Description

### SHUT BAG HOLDER - 70 l

Trolley suitable for carrying one 70 l bag (not included).  
Stainless steel AISI 304 structure.  
Equipped with swivelling castors with brake and bumpers.  
For use with standard 70 l plastic bag.  
Can use also wrapping bag (27462, 27466, 45908)

Made in Italy.