

CLOSED LAUNDRY TROLLEY - painted steel



| | |
|---------------|-------------------|
| Product Code | 45925 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | No |
| CH-REP | No |
| EAN/UPC | 8023279459258 |

Description

CLOSED LAUNDRY TROLLEY - painted steel

Trolleys for collection of dirty laundry and distribution of clean linen to the patient.

Stainless steel or painted steel frame.

Equipped with 4 Ø 125 mm anti-thread swiveling castors (2 with brake), 2 rectangular supports with rubber bag stop ring, 2 lids and 3 shelves of which 2 adjustable.

Suitable for 70 litres waste bags.

Size: 1330x650xh 1200 mm (45925) 1330x650xh 1140 mm (45926)

Made in Italy.