

CLEAN LINEN TROLLEY 2 doors + 1 drawer - small



|               |                   |
|---------------|-------------------|
| Product Code  | 45931             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | No            |
| CH-REP  | No            |
| EAN/UPC | 8023279459319 |

Description

CLEAN LINEN TROLLEYS 2 DOORS

Clean linen trolley - small  
Size (mm): 670x640xh 1280  
Drawer size (mm): 600x480xh 150

Solid and sturdy trolleys for distribution of clean linen and patient hygiene.  
Steel painted with antimicrobial epoxy powders to obtain maximum lightness.  
Upper part has large support surface in technopolymer plastic material.  
Integrated ergonomic push handle.

Supplied with 2 doors with 2 internal shelves and one drawer made of antimicrobial epoxy powder coated steel. Key closures. 4 castors Ø 125 mm, anti-wire, pivoting, two of which with lever locking system and reinforced plastic support.

Made in Italy.